

Aespire* View

Clarity of view...Clarity of decision

Features

- Large 12-inch color ventilator display with color waveform and alarm message indicators
- Lightweight and compact for easy maneuverability
- Optional Total Fresh Gas modules for electronic measurement and display of the total fresh gas flow and individual gas flow (O₂, Air, N₂O)
- Optional integrated auxiliary O₂ flowmeter and suction control
- Enhanced monitor integration capabilities with our GE Healthcare Anesthesia Monitor and Compact Monitor

Advanced Ventilation: 7900 SmartVent*

- Maximum versatility for full patient range
- Ventilation Modes:
 - VCV (Volume Control)
 - PCV Pressure Control (optional)
 - PSVPro* (Pressure Support with Apnea Backup) (optional)
 - SIMV/PSV (Synchronized Intermittent Mandatory Ventilation with pressure support) (optional)
 - SIMV-PC (Synchronized intermittent Mandatory Ventilation with pressure control) (optional)
 - PCV-VG (Pressure Controlled Ventilation – Volume Guaranteed) (optional)
- Electronic PEEP
- Automatic fresh gas flow (tidal volume) compensation
- Cardiac bypass case mode
- Direct access to ventilator parameter settings



Aespire* View shown with GE Tec 7* Vaporizer and CARESCAPE Monitor B650

- Pressure waveform for visual reference on a breath-by-breath basis
- Smart alarms direct user to specific problems and affected parameters
- Inspired oxygen monitoring

Advanced Breathing System (ABS)

- Easy to clean, autoclavable, latex-free
- Fast response – exceptional for low flow anesthesia
- Easy removal – no tools required
- Integrated design – less parts and connections helps reduce potential for leaks and misconnects
- One step bag/vent switch turns ventilator on/off
- Optional CO₂ bypass with electronic detection of “absorber off”
- Optional EZChange and condensor



Physical Specifications

Dimensions

Height:	136 cm/53.5 in
Width:	75 cm/29.5 in
Depth:	74 cm/29.1 in
Weight:	approximately 136 kg/300 lbs

Top shelf

Weight limit:	34 kg/75 lbs
Width:	65 cm/26 in
Depth:	43 cm/15.75 in

Work surface

Height:	83.8 cm/33 in
Size:	2160 cm²/334 in²

Folding side shelf (optional)

Height:	87.5 cm/34.5 in
Width:	26.5 cm/10.4 in
Depth:	31.5 cm/12.4 in
Weight limit:	11.3 kg/25 lbs

DIN rail

Side of machine: 34.5 cm/13.6 in

Drawers (internal dimensions)

Height:	17.5 cm/6.9 in
Width:	33 cm/13 in
Depth:	26.5 cm/10.4 in

Absorber bag arm (optional)

Arm length:	30.5 cm/12 in
Bag arm height (adjustable):	87 cm/34.3 in 104 cm/40.9 in

Casters


Diameter:	12.5 cm/5 in
Brakes:	Individual locking

Ventilator screen

Display size	31 cm/12.1 in (diagonal)
--------------	--------------------------



Ventilator Operating Specifications



Ventilation operating modes

VCV
PCV
SIMV/PSV
PSVPro* with Apnea backup
SIMV-PC
PCV-VG

Ventilator (VT) parameter ranges

Tidal volume range: 20 to 1500 mL (Volume Control SIMV/PSV and PCV-VG modes)
5 to 1500 mL (Pressure Control Mode)

Incremental settings: 20 to 50 mL (increments of 1 mL)
50 to 100 mL (increments of 5 mL)
100 to 300 mL (increments of 10 mL)
300 to 1000 mL (increments of 25 mL)
1000 to 1500 mL (increments of 50 mL)

Minute volume range: 0 to 99.9 L/min

Pressure (P_{inspired}) range: 5 to 60 cm H₂O (increments of 1 cm H₂O)

Pressure (P_{limit}) range: 12 to 100 cm H₂O (increments of 1 cm H₂O)

Pressure (P_{support}) range: Off, 2 to 40 cm H₂O (increments of 1 cm H₂O)

Rate: 4 to 100 breaths per minute for PCV-VG, Volume Control and Pressure Control vent modes
2 to 60 breaths per minute for SIMV/PSV, PSVPro and SIMV vent modes (increments of 1 breath per minute)

Inspiratory/expiratory ratio: 2:1 to 1:8 (increments of 0.5)

Inspiratory time: 0.2 to 5.0 seconds (increments of 0.1 seconds) (SIMV and PSVPro vent modes)

Trigger window: 0 to 80% (increments of 5%)

Flow trigger: 0.2 to 1.0 L/min (increments of 0.2 L/min)
1 to 10 L/min (increments of 0.5 L/min)

Inspiration termination level: 5 to 75% (increments of 5%)

Backup mode delay: 10 to 30 seconds (increments of 5 seconds)

Positive End Expiratory Pressure (PEEP)

Type: Integrated, electronically controlled

Range: OFF, 4 to 30 cm H₂O (increments of 1 cm H₂O)

Ventilator performance

Pressure range at inlet: 240 kPa to 700 kPa/
35 psig to 100 psig

Peak gas flow: 120 L/min + fresh gas flow

Flow valve range: 1 to 120 L/min

Fresh gas flow compensation

Ventilator monitoring

Expiratory minute volume range: 0 to 99.9/min

Expiratory tidal volume range: 0 to ≥ 1500 mL

O₂ %: ≤ 5 to 110%

Peak pressure: -20 to 120 cm H₂O

Mean pressure: -20 to 120 cm H₂O

Plateau pressure: 0 to 120 cm H₂O

Pressure waveform sweep speed: 4 to 25 breaths per minute (0 to 16 seconds)
26 to 75 breaths per minute (0 to 8 seconds)
75 breaths per minute (0 to 4 seconds)

Ventilator accuracy

Delivery/monitoring accuracy

Volume delivery:	> 210 mL = better than 7% < 210 mL = better than 15 mL < 60 mL = better than 10 mL
Pressure delivery:	$\pm 10\%$ or ± 3 cm H ₂ O
PEEP delivery:	± 1.5 cm H ₂ O
Volume monitoring:	> 210 mL = better than 9% < 210 mL = better than 18 mL < 60 mL = better than 10 mL
Pressure monitoring:	$\pm 5\%$ or ± 2 cm H ₂ O

Alarm settings

Tidal volume (TV _{exp}):	Low: OFF, 5 to 1500 mL High: 20 to 1600 mL, OFF
Minute volume (Mv _{exp}):	Low: OFF, 0.1 to 10 L/min High: 0 to 30 L/min, OFF
Inspired oxygen (FiO ₂):	Low: 18 to 99% High: 21 to 99%, OFF
Apnea alarm:	Mechanical ventilation ON: < 5 mL breath measured in 30 seconds Mechanical ventilation OFF: < 5 mL breath measured in 30 seconds
Low airway pressure:	4 cm H ₂ O above PEEP
High pressure:	12 to 100 cm H ₂ O (increments of 1 cm H ₂ O)
Sustained airway pressure:	<i>Mechanical ventilation ON:</i> $P_{\text{limit}} < 30$ cm H ₂ O, the sustained limit is 6 cm H ₂ O $P_{\text{limit}} 30$ to 60 cm H ₂ O, the sustained limit is 20% of P_{limit} $P_{\text{limit}} > 60$ cm H ₂ O, the sustained limit is 12 cm H ₂ O <i>PEEP and mechanical ventilation ON:</i> Sustained limit increases by PEEP minus 2 cm H ₂ O <i>Mechanical ventilation OFF:</i> $P_{\text{limit}} \leq 60$ cm H ₂ O, the sustained limit is 50% of P_{limit} $P_{\text{limit}} > 60$ cm H ₂ O, the sustained limit is 30 cm H ₂ O
Subatmospheric pressure:	Paw < -10 cm H ₂ O
Alarm silence countdown timer:	120 to 0 seconds

Ventilator components

Flow transducer

Type:	Variable orifice flow sensor
Dimensions:	22 mm OD and 15 mm ID
Location:	Inspiratory outlet and expiratory outlet
(optional autoclavable sensor available)	

Oxygen Sensor

Type:	Galvanic fuel cell
Life Cycle	Approximately 12 months (Dependent on usage)

Anesthetic agent delivery

Delivery

Vaporizers:	Tec 6 Plus, Tec 7
Number of positions:	2
Mounting:	Tool-free installation Selectatec® manifold interlocks and isolates vaporizers



Tec 6 Plus Vaporizer



Tec 7 Vaporizer

Electrical specifications

Current leakage

100/120 V:	< 300µA
220/240 V:	< 500µA

Power and battery backup

Power input:	100-120 Vac, 50/60 Hz 220-240 Vac, 50/60 Hz
Backup power:	Demonstrated battery backup time under typical operating conditions is 90 minutes when fully charged
Battery type:	Internal rechargeable sealed lead acid
Power cord:	Length: 5 m/16.4 ft or 3m /9.8ft Rating: 10A @ 220 Vac or 15A @ 120 Vac

Communication port

Serial interface:	RS-232 compatible port
-------------------	------------------------

Inlet/outlet modules

	120 V
System circuit breakers:	15A
Outlets:	4 outlets on back, 3-2A, 1-3A individual breakers with isolation transformer

Pneumatic specifications

Auxiliary common gas outlet

Connector:	ISO 22 mm OD and 15 mm ID
------------	---------------------------

Gas supply

Pipeline input range:	280 kPa to 600 kPa/41 psi to 87 psi
Pipeline connections:	DISS-male, DISS-female, AS4059, FS90-116, SS 87 524 30, or NIST (ISO 5359). All fittings available for O ₂ , N ₂ O, and Air, and contain pipeline filter and check valve.
Cylinder input:	Pin indexed in accordance with CGA-V-1 or DIN (nut and gland); contains input filter and check valve Note: Maximum 3 cylinders; two inboard mounted, one outboard mounted.

Primary regulator diaphragm minimum burst pressure:	2758 kPa/400 psig
Primary regulator nominal output:	≤338 kPa/49 psig Pin indexed cylinder connections ≤407 kPa/59 psig DIN cylinder connections

O₂ controls

Method:	Proportionate decrease of N ₂ O with reduction in O ₂ pressure
Supply failure alarm:	Range: 193 kPa to 221 kPa/28 psig to 32 psig Sounds at maximum volume every 10 seconds
O ₂ flush:	Range: 35 to 75 L/min

Flowmeters

O ₂ ranges:	0.05 to 0.95 L/min and 1.0 to 15.0 L/min; Minimum O ₂ flow: 50 mL/min ±25 mL	
N ₂ O ranges:	0 to 0.95 L/min and 1.0 to 10.0 L/min	
Air range:	0 to 0.95 and 1.0 to 15.0 L/min	
Calibration:	Percent of full scale flow	Accuracy (% of flowrate)
	100	±2.5%
	90	±2.5%
	80	±2.6%
	70	±2.7%
	60	±2.9%
	50	±3.1%
	40	±3.4%
	30	±4.0%
	20	±5.0%
	10	±8.1%

Calibration conditions:* 20°C/68°F, 101.3 kPa/760 mmHg

* Different breathing circuit pressures, barometric pressures or temperatures change flowtube accuracy.

Total flow sensor

Accuracy: Greater of ± 25 sccm (smlpm) or $\pm 6.0\%$ of measured value

Calibration conditions: $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$, 69.7 ± 2 psi

Note: Different breathing circuit pressures, barometric pressures or temperature changes total flow sensor accuracy

Hypoxic guard system

Type: Mechanical Link-25

Range: Provides a nominal minimum 25% concentration of oxygen in $\text{O}_2/\text{N}_2\text{O}$ mixture

Environmental specifications

System operation

Temperature: 10° to 40°C / 50° to 104°F

Humidity: 15 to 95% relative humidity

Altitude: -440 to 3565 m/ 500 to 800 mmHg

System storage

Temperature: -25° to 65°C / -13° to 149°F

Humidity: 15 to 95% relative humidity

Altitude: -440 to 5860 m/ 375 to 800 mmHg

Oxygen cell storage: -15° to 50°C / 5° to 122°F
15 to 95% relative humidity
 500 to 800 mmHg

Electromagnetic compatibility

Immunity: Complies with all requirements of EN 60601-1-2

Emissions: CISPR 11 group 1 class B

Approvals: UL 60601-1,
CSA C22.2 #601.1
EN/IEC 60601-1
CE

Breathing circuit specifications

Operational modes

Breathing circuit is circle mode only

Carbon dioxide absorbent canister

Absorbent capacity: 800 g

Integrated expiratory limb water reservoir

Ports and connectors

Exhalation: 22 mm OD ISO 15 mm ID taper

Inhalation: 22 mm OD ISO 15 mm ID taper

Bag port: 22 mm OD

Pressure gauge

Scale range: -2 to 10 kPa/ -20 to 100 cm H₂O

Bag-to-Ventilator switch

Type: Bi-stable

Control: Controls ventilator and direction of breathing gas within the circuit

Integrated Adjustable Pressure Limiting (APL) valve

Range: 0.5 to 70 cm H₂O

Tactile knob indication at: 30 cm H₂O and above

Adjustment range of rotation: 0.5 to 30 cm H₂O (0 to 230°)
 30 to 70 cm H₂O (230 to 330°)

Materials

All materials in contact with exhaled patient gases are autoclavable, except disposable flow sensors and O_2 cell. (Autoclavable flow sensors optional).

All materials in contact with patient gas are free of natural rubber latex.

Breathing circuit parameters

Compliance: Bag mode: 1.82 mL/cm H₂O

Mechanical mode: Automatically compensates for compression losses within the absorber and bellows assembly

Circuit volume: 2.7 L Vent Mode (including absorber)
 1.2 L Bag Mode

Expiratory resistance:	Flow rate	Pexp Bag Mode Pressure drop	Pexp Vent Mode Pressure drop
	10 L/min	0.78 cm H ₂ O	0.77 cm H ₂ O
	30 L/min	1.59 cm H ₂ O	1.71 cm H ₂ O
	60 L/min	3.48 cm H ₂ O	3.88 cm H ₂ O

Note: With patient circuit and wye piece add $+0.89$ cm H₂O

Anesthetic gas scavenging

System Type	Outlet Connector	Hospital waste gas disposal system requirements
Active Adjustable Flow, High Vacuum	DISS EVAC	12inHg (305 mmHg) minimum vacuum @ 36 LPM flow Extract Flow limited to 30 LPM @ 12 inHg vacuum
Active Low Flow, High Vacuum	DISS EVAC	12 inHg (305 mmHg) minimum vacuum @ 36 LPM flow
Active Low Flow, Low Vacuum	12.7 mm Barb	36 LPM minimum flow

© 2011 General Electric Company – All Rights Reserved.

General Electric Company reserves the right to make changes in specification and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

*Aespire, Tec, and PSVPro are trademarks of General Electric Company.

Datex-Ohmeda, Inc., a General Electric company, doing business as GE Healthcare.

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our “healthymagination” vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality around the world.

Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

GE Healthcare
3030 Ohmeda Drive
P.O. Box 7550
Madison, WI 53707-7550
U.S.A.

www.gehealthcare.com



imagination at work