

For Veterinary use only



All-in-one Mobile C-arm with FPD for Veterinary Use

VETERA ACE





Key Features

Experience Innovation in Space, Imaging, and Procedures with the Integrated FPD

All-in-One Body

Compact Design for Optimal Use of Limited Space.

Free Space 860mm

Enhanced Free Space for Unrestricted Procedures.

High-Speed Image Processing

GENORAY's advanced low-latency technology minimizes image delay, delivering fast and accurate images that maximize procedural efficiency.

Low-Dose Technology for Animals and Users

GENORAY's low-dose technology reduces radiation while maintaining high image quality, protecting animals and minimizing long-term exposure risks for users.

Product Overview



All-In-One Body Design



3.2kW HFG



Free Space 860mm



20x20cm CMOS
/ 21x21cm TFT




27" FHD
Touch Monitor



Digital Subtraction
Angiography


The All-in-One body, which integrates the main unit and key components, operates without the need for an additional cart, optimizing space in the procedure room.

VETERA ACE



DNR (Dynamic Noise Reduction)


Dynamic Noise Reduction eliminates image noise, delivering sharper and more reliable images.



Dose Level Adjustment

Optimized dose adjustments ensures safe use in various dose modes, while consistently delivering clear images for various animal sizes and procedures.

Comprehensive Veterinary Imaging



Style Template Function


Offers three customizable imaging styles, allowing users to choose their preference and create an optimal imaging environment.



Snap Shot

While in real-time fluoroscopy, users can capture their desired still image in high resolution.

Versatile Imaging Across Species and Sizes



Real-Time Image Comparison – Swap

Freely switch between live and reference images, to enable an accurate comparative diagnosis and localization for procedures.




Image Scroll Function

Stored images can be reviewed quickly while scrolling, making it much more efficient.


Contrast-Enhanced Fluoroscopic Imaging

* DSA Option Available



Dimensions

* Minimum Installation Space



GENORAY

GENORAY Co.,Ltd.

T +82-31-5178-5500
E genoray@genoray.com
W www.genoray.com

GENORAY EU GmbH

T +49-30-50969498
E info@genorayeuropa.com
W europe.genoray.com

GENORAY Türkiye

T +90-216-463-57-37, +90-216-463-57-53
E genorayturkiye@genoray.com
W turkiye.genoray.com

GENORAY America Inc.

T +1-855-436-6729
E inquiry@genorayamerica.com
W www.genorayamerica.com

GENORAY Japan

T +81-45-620-4971
E Info@genorayjapan.co.jp
W www.genorayjapan.co.jp

GENORAY Shanghai

T +86-021-6785-8538
E genoraychina@genoray.com
W www.genoray.com.cn



* This X-ray unit may be dangerous to operator unless safe exposure factors and operating instructions are observed.

© GENORAY CO.,LTD. Ver. VETERA_ACE-Rev_00_260126_ENG