



VistOR PRO, PR, EX Surgical Lights



The guiding lights for all patient care

Whether you're performing surgery, procedures or exams, it's absolutely critical that the lighting conditions be optimal so you can visualize your work with complete accuracy. Nothing less is acceptable. At Nuvo, our sole focus is making that happen.

TheVistOR lights are significantly brighter, whiter and cooler than any previous technology. So you get the most vivid, accurate colors with remarkable shadow control — and virtually no heat emission — which make our lights more comfortable to work under than any other lights. All at a price you won't sweat over.

VistOR At A Glance

- 4,300° K pure white illumination
- 100,000/75,000/55,000 lux output
- 50,000 hour-rated LED life
- CRI (Color Rendering Index) of 97
- Superior shadow control
- Near-zero radiated heat
- Wide array of mounts for ultimate flexibility



Getting a better handle on things.

Controlling the brightness and position of your lights is critical. With the VistOR LED Series, it's incredibly simple. Just reach up and press the center button on our sterilizable multi-stage dimming control handle to adjust brightness.

Perfect balance and greater shadow control.

The arm of each light is precision-balanced and articulated to allow the greatest reach and range of positioning with surprising ease. Multiple lighting elements and a large diameter platform, especially on the VistOR Pro, give you an amazing amount of light output along with remarkable shadow control, so you can visualize the site under the most optimal conditions possible.



multi-stage dimming and illumination field are easily controlled with the sterilizable handle.



Engineering meets Science meets Surgery.

No one else in the industry has gone to the extremes of engineering as we have. But then, bringing perfection to medial lighting is what we do. In developing our superb VistOR Series, we sought to create the most pure, uniform white light possible, taking into account not only how LEDs operate but also how the human eye perceives color. The result is a remarkable optics technology that sets a new standard in surgical illumination.



To compensate for inherent LED color imperfections, we developed an ingenious lens that merges all emitted colors into one focused beam of very white light, affording you unparalleled color accuracy.



The VistOR LEDs incorporate a special front diffuser lens system that individually re-directs each beam into a single, coherent, even light field.

VistOR Pro Series

Three incomparable lights with limitless lighting possibilities.

Our VistOR Series offers a selection of three different-sized lights. Choose the best one to fit your operating, procedure or exam rooms.

VistOR Pro – Adding unsurpassed brilliance in the OR.

Ideal within the surgical environment, the VistOR Pro offers unsurpassed longevity and flexibility in a package that's incredibly affordable. The large diameter LED array provides excellent shadow control – even with just a single light head. Pressing the center button on the sterilizable handle allows 5-stage dimming plus ambient light mode (for minimally-invasive procedures) for easy management of brightness within the illumination field.

- 100,000 lux at 1 meter
- 4,300° K pure white illumination



VistOR Pro
Single Mount



VistOR PR – The best choice for brighter procedure outcomes.

The ideal and surprisingly affordable, mid- sized light for minor surgery, procedures and exams. The VistOR PR delivers maximum brightness, sharpness and color accuracy. So you can visualize your work in the most precise way. 3-stage dimming via the sterilizable handle lets you easily control brightness within the sterile field.

- 75,000 lux at 1 meter
- 4,300° K pure white illumination



VistOR PR
Dual Mounts




With Sterilizable Handle
and Control Grip

VistOR EX – LED Excellence for exams and diagnostics.

The VistOR EX is the perfect fit and offers the greatest value in an examination and diagnostic light. It provides maximum brightness, sharpness and color accuracy. So you can visualize your work in the most precise way.

- 55,000 lux at 24" (60 cm)
- 4,300° K pure white illumination





Dual Ceiling Mount

Single Ceiling Mount

Wall Mount

Portable

Design your lighting to fit your exact needs.

The VistOR Series gives you total flexibility in how you design the lighting in your OR, ER, ICU, clinic or office. Featuring a longer, completely reengineered arm (now 66"), the VistOR Pro, PR and EX are available in single ceiling, dual ceiling, wall, and portable mounting options.

Optional 7-hr Battery Back-Up for portable and ceiling/wall mounts (VistOR PR and VistOR Pro)



NUVO - Operating in the patient's best interest

For nearly 20 years, our sole focus at Nuvo has been to provide the medical field with the highest quality, most cost-effective surgical lights, tables, organizing systems and other products for hospital operating rooms. All of which can greatly help you make a difference in the care of your patients.

Check the specs at www.nuvosurgical.com or contact your distributor.



(800) 663-1152

www.nuvosurgical.com



Proudly made in the USA



Nuvo is a division of Medical Illumination
19749 Dearborn Street, Chatsworth, CA

