

*Redesigned*

# General Transport Stretcher PT1001

**750**

POUND

Weight  
Capacity



## STANDARD FEATURES

- 750 lb. weight capacity
- Height adjustment range: 22"-35" (bottom of mattress)
- Ergonomically located oxygen tank holder
- Effortless steering, pivoting, and cornering with 8" Tente casters (1 directional & 3 swiveling)
- 6 brake/steer pedals
- Dual pedestal lift mechanism
- Lift/Trendelenburg pedals on both sides
- Trendelenburg & reverse Trendelenburg +/- 18°
- Dual Fowler assisted lift back rest with 90° adjustability
- Fold forward side rails with built in push handles
- 29" or 31" bed width
- 4" comfort mattress
- 2 hook removable telescopic IV pole
- 6 IV pole sockets available
- Large storage area underneath for patient belongings
- Full perimeter bumper to protect your walls and fixtures
- 4 restraint belt locations

## OPTIONAL FEATURES

- 5th wheel add-on kit
- Fold down push handles
- Knee gatch
- Paper roll dispenser kit \*
- Detachable foot board \*
- Envelope style rail bumper pads \*

\* Additional lead time required (7-10 business days)

## PRODUCT SPECIFICATIONS

### PT1001

- 29" bed width
- Overall dimensions: 83 5/16" long x 33 9/16" wide
- Patient surface dimensions: 75 7/16" long x 29 1/8" wide

### PT1001-31

- 31" bed width
- Overall dimensions: 83 5/16" long x 35 9/16" wide
- Patient surface dimensions: 75 7/16" long x 31 1/8" wide

MAC Medical, Inc. is pleased to announce the completion of our redesigned General Transport Stretcher. The PT1001 allows for easy and comfortable patient transfer with an increased weight capacity of up to 750 lbs. Economical and versatile, the PT1001 is an ideal solution for any department where cost efficient quality is required.


For questions or to place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at [sales@macmedical.com](mailto:sales@macmedical.com).

**MADE** IN THE **USA**

**mac**  
**MEDICAL**

[www.macmedical.com](http://www.macmedical.com)

## Features




Effortless steering, pivoting, and cornering with 8" Tente casters (1 directional & 3 swiveling)

Dual Fowler assisted lift back rest with 90° adjustability




Ergonomically located oxygen tank holder




Large storage area for patient belongings




6 brake/steer pedals (two on each side & one on each head and foot end)



5th wheel add-on kit (optional)




Trendelenburg & reverse Trendelenburg +/- 18°




Height adjustment range: 22"-35" (bottom of mattress)

## Dimensions


### PT1001



### PT1001-31



### Height Adjustment Ranges (Both Models)



**Customer needs are our first priority.**

MAC Medical, Inc.  
Printed in USA  
Publication No. FFU-16 Rev B  
October 2016

Information regarding this product is subject to change without prior notice.



**Corporate Office**  
325 West Main Street  
Belleville, IL 62220

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • Fax: 618-476-3337  
Email: sales@macmedical.com • Website: www.macmedical.com