


New

5th Wheel Transport Stretcher PTF1000



750
POUND
Weight
Capacity

Retractable
Fifth Wheel

MAC Medical, Inc. is pleased to announce the arrival of the new PTF1000 5th Wheel Transport Stretcher. The retractable center fifth wheel combines effortless directional movement and cornering with 4 swiveling wheels for ease of use, control, and safety.

For questions or to place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at sales@macmedical.com.

**MADE
IN THE USA**

Customer needs are our first priority.

STANDARD FEATURES

- Retractable center fifth wheel provides effortless directional movement and cornering
- 750 lb. weight capacity
- Height adjustment range: 22"-35" (bottom of mattress)
- Ergonomically located oxygen tank holder
- Effortless steering, pivoting, and cornering with 8" Tente casters (4 swiveling)
- 6 brake/steer pedals
- Dual pedestal lift mechanism
- Lift/Trendelenburg pedals on both sides
- Trendelenburg & reverse Trendelenburg +/- 18°
- Dual Fowler assisted lift back rest with 90° adjustability
- Fold forward side rails with built in push handles
- 29" or 31" bed width
- 4" comfort mattress
- 2 hook removable telescopic IV pole
- 6 IV pole sockets available
- Large storage area underneath for patient belongings
- Full perimeter bumper to protect your walls and fixtures
- 4 restraint belt locations

OPTIONAL FEATURES

- Fold down push handles
- Knee gatch
- Paper roll dispenser kit *
- Detachable foot board *
- Envelope style rail bumper pads *

* Additional lead time required (7-10 business days)

PRODUCT SPECIFICATIONS

PTF1000

- 29" bed width
- Overall dimensions: 83 5/16" long x 33 9/16" wide
- Patient surface dimensions: 75 7/16" long x 29 1/8" wide

PTF1001-31

- 31" bed width
- Overall dimensions: 83 5/16" long x 35 9/16" wide
- Patient surface dimensions: 75 7/16" long x 31 1/8" wide

mac
MEDICAL

www.macmedical.com

Features




Retractable center fifth wheel provides effortless directional movement and cornering

Dual Fowler assisted lift back rest with 90° adjustability



Ergonomically located oxygen tank holder




Large storage area for patient belongings




6 brake/steer pedals (two on each side & one on each head and foot end)



Effortless steering, pivoting, and cornering with 8" Tente casters (4 swiveling)




Trendelenburg & reverse Trendelenburg +/- 18°




Height adjustment range: 22"-35" (bottom of mattress)

Dimensions


PTF1000



PTF1000-31



Height Adjustment Ranges (Both Models)



Customer needs are our first priority.

MAC Medical, Inc.
Printed in USA
Publication No. FFU-01 Rev B
October 2016

Information regarding this product is subject to change without prior notice.



Corporate Office
325 West Main Street
Belleville, IL 62220

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • Fax: 618-476-3337
Email: sales@macmedical.com • Website: www.macmedical.com