

# Prep & Pack Workstations

Ergonomical and versatile workstation for all body types

Ergonomically  
Friendly

Designed with ergonomics in mind, the Prep & Pack Workstations feature a modular design, allowing for multiple configurations and adjustability for all body types. Choose from the many optional features including adjustable height (manual or electronic), overhead light, articulating magnification light, CPU holder, electrical strip, and more.

## STANDARD FEATURES

- 100% stainless steel construction
- Modular design allows for multiple configurations
- Integrated footrest
- 4" casters (2 swivel & 2 locking)

## OPTIONAL FEATURES

- **Ergonomically friendly** adjustable heights (manual or electronic)
  - Adjustable height range: 36"–44"
- Bin rail
- Overhead light
- Articulating magnification light
- Electrical strip
- CPU holder
- Pull out keyboard tray
- Magnetic dry erase board
- Adjustable flat shelf for riser frame
- Adjustable angled shelf for riser frame
- Lower shelf
- Side mount waste basket holder with waste basket
- Spool holder
- Drawer boxes
  - Single 6" drawer box
  - Double 6" drawer box



**TWS-3672**  
(shown with multiple optional features)

## Standard Models

Base Model #	Overall Dimensions (height at lowest setting)
TWS-3060	60"W x 32.25"D x 81"H
TWS-3072	72"W x 32.25"D x 81"H
TWS-3660	60"W x 38.25"D x 81"H
TWS-3672	72"W x 38.25"D x 81"H

*Note: Additional sizes available for specific requirements. Please contact us for further assistance.*

*Note: Model # will change based on features added.*

**MADE IN THE USA**

**mac**  
**MEDICAL**

**Corporate Office**  
325 West Main Street  
Belleville, IL 62220

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550  
Toll Free: 877-828-9975  
Fax: 618-476-3337  
sales@macmedical.com  
www.macmedical.com

MAC Medical, Inc.  
Publication No. FFU-06 Rev D  
August 2017  
Information regarding this product is subject to change without prior notice.