

# PG SERIES

## *Exam Table*

- Designed for fast paced, diverse treatments
- Open base design for great ergonomics
- Multi-purpose PG350 converts from a table to a chair
- ADA compliant height range



PG350  
shown in Heron fabric



PG200



PG300



PG350 shown with upgrades



ISO 13485:2003  
FM 584901

made in the USA  
with US & imported parts



**OAKWORKS<sup>®</sup>**  
**MEDICAL**  
*Position for Success*

# PG SERIES

## Exam Table

The PG Series Exam Tables are designed for busy environments where they may need to be moved often to permit large pieces of equipment to occupy the same space. Casters come standard on this model, and several tops are available along with many options and accessories to support multiple procedures.



ZERO|GRAVITY  
3 Section Top

### MODEL PG350

#71953-T\* 30" W 3" Comfort  
#72422-T\* 30" W 4" Comfort

|          |                                    |
|----------|------------------------------------|
| Motion 1 | Electric Height: 22"-38" (56-97cm) |
| Motion 2 | Electric Tilt: 15°/15°             |
| Motion 3 | Electric Backrest: 0-75°           |
| Motion 4 | Manual 4 Position Leg rest: 0-90°  |



ZERO|GRAVITY  
Backrest Top

### MODEL PG300

#71952-T\* 30" W 3" Comfort  
#72420-T\* 30" W 4" Comfort

|          |                                    |
|----------|------------------------------------|
| Motion 1 | Electric Height: 22"-38" (56-97cm) |
| Motion 2 | Electric Tilt: 15°/15°             |
| Motion 3 | Electric Backrest: 0-75°           |



Flat Top

### MODEL PG200

#71948-T\* 30" W 3" Comfort  
#72419-T\* 30" W 4" Comfort

|          |                                    |
|----------|------------------------------------|
| Motion 1 | Electric Height: 22"-38" (56-97cm) |
| Motion 2 | Electric Tilt: 15°/15°             |

| FEATURES          |  |
|-------------------|--|
| Top               | 30" (76cm) wide x 72" (183cm) long                     |
| Foot Control      | Controls all powered movements                         |
| Weight Capacity   | 550 lb (250kg)   |
| Product Weight    | 159 lb (163kg)   |
| Padding           | 3" (8cm) or 4" (10cm) Comfort Foam™                    |
| Fabric            | Premium, PVC free TerraTouch™                          |
| Casters           | Heavy-duty 4" Locking                                  |
| Paper Roll Holder | Accommodates 21" (53cm) roll (roll not included)       |
| Warranty          | 3 years parts, 1 year labor                            |
| Safety Listings   | FDA listed, CE marked, cETLus approved, CB Certificate |

| OPTIONS <i>must be ordered with table.</i> |   | #          |
|--|---|------------|
| Accessory Mounting Plate                   | Allows the use of accessories with your table                                     | 62123      |
| Central Base Locking System                | Not available with Patient Lift Compatible Base                                   | 61777      |
| Non-USA electric system                    | please contact us for order details   | 73134      |
| Patient Lift Compatible Base               | Not available with Central Base Locking System                                    | 65893      |
| Procedure Tray                             | 12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters) (PG350 30"W)         | 73042      |
| Side Rails (set of 2)                      | 23.5" (60cm) long x 6" (15cm) tall (from surface of table)                        | 61781      |
| Stirrups (Retractable Foot Rests)          | Pull out from under seat section, fully adjustable (PG350)                        | 71001      |
| ACCESSORIES <i>order any time.</i>         |   | #          |
| T-rail Adapter                             | Attaches to the backrest for use with USA format T-rail attachable accessories.   | 61902      |
| Detachable Arm board                       | Use with T-Rail Adapter - Rotates over 140°<br>5" (13cm) wide x 20.5" (52cm) long | 62192      |
| Universal Table Extender with 3" Padding   | 21" (53cm) W x 12.5" (32cm) L (requires mounting plate)                           | 68697      |
| Universal Table Extender with 4" Padding   | 21" (53cm) W x 12.5" (32cm) L (requires mounting plate)                           | 68735      |
| QuickLock™ Head Rest (black) & AeroCel Pad | Comfortable, adjustable prone or supine support. (requires mounting plate)        | PKG4186-T* |

\*Fabric color code must be specified.

## CENTRAL BASE LOCKING SYSTEM #61777

Easily secure your table into place with the central base locking system.



OPTIONAL

## SIDE RAILS #61781

Increase stability and security for your clients with a pair of side rails. Fold down out of the way when not in use. 23.5" (60cm) long x 6" (15cm) tall (from surface of table)



OPTIONAL

## STIRRUPS #73175

Adjust to a wide range of angles & lengths to accommodate any patient or procedure. Retractable design stores inside the table.



(PG350 only)  
OPTIONAL

## HAND CONTROL #enquire for part number

Add a hand control to control all functions.



OPTIONAL

## PROCEDURE TRAY #73042

Retractable, removable stainless steel. 12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters)



OPTIONAL

## PATIENT LIFT #65893 COMPATIBLE BASE

This option allows the use of patient lifts and easy wheel chair transfers to minimize or eliminate caregiver injuries associated with patient transfers. Must be ordered with table.



OPTIONAL

## T-RAIL ADAPTER #61902

Multiple locations allow you to choose the best spot for this easily installed & removed accessory capable of handling up to 150 lb (68kg). Can be used with any American standard accessory. 8" (20cm) long



ACCESSORY

## QUICKLOCK™ HEAD REST & AERO·CEL PAD #PKG4186-T\*

The platform allows precise cervical flexion and height positioning over a wide range for both prone and supine treatments.



ACCESSORY

## UNIVERSAL TABLE EXTENDER #68697-T\* 3" padding #68735-T\* 4" padding


Perfect as a head rest or foot rest, it has two rods that fit into outlets on the table. Extends table to 84". 21" (53cm) W x 12.5" (32cm) L



ACCESSORY

### TERRATOUGH<sup>™</sup> FABRIC

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables and passes testing for Cytotoxicity, Skin Irritation and Sensitization as well as CA TB117 Fire Retardance testing. For disinfecting purposes you may use Protex, MadaCide, MadaWipe, Accell TB, Virox<sup>®</sup> (Canada), PDI Duper Sani-Cloth Germicidal wipe or a diluted 1/10 bleach solution.



*We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.*

### PADDING



**3" Comfort Foam<sup>™</sup>**  
*More comfort and support*



**4" Comfort Foam<sup>™</sup>**  
*Superior comfort and support*

CONTACT:



*View our full  
line of Medical  
products #72703*



Scan this code  
with your  
smart phone