

Compare the **Key Features** of our Power Exam Tables

Key Features

Key Features	proglide300+	proglide300				
	4040-650-300	4070-650-300	4040-650-200	4070-650-200	4040-650-100	4070-650-100
Mobility & Efficiencies						
OneTouch Wheelbase® System	✓	✓	-	-	-	-
Hand Control with OneTouch Patient Positioning® System 2 Programmable Positions & Return-to-Home	✓	✓	✓	✓	-	-
Two function foot control with Hi-Lo and power back	✓	-	✓	-	✓	-
One function foot control with Hi-Lo	-	✓	-	✓	-	✓
Power Hi-Lo + Power Backrest	✓	-	✓	-	✓	-
Power Hi-Lo with Air-Spring Assisted Manual Backrest	-	✓	-	✓	-	✓
Dual electro-mechanical actuator system for efficient speed	✓	✓	✓	✓	✓	✓

Accessibility	proglide300+	proglide300				
	4040-650-300	4070-650-300	4040-650-200	4070-650-200	4040-650-100	4070-650-100
650lb (294.84kg) Patient Weight Capacity	✓	✓	✓	✓	✓	✓
Low-access ADA Height 18.5" (±0.5") – 46.9cm (±1.27cm)	✓	✓	✓	✓	✓	✓
Pelvic tilt (7°)	✓	✓	✓	✓	✓	✓
Drain Pan: Stainless steel 1.5" deep (3.81cm)	✓	✓	✓	✓	✓	✓
Large front storage drawer	✓	✓	✓	✓	✓	✓
Side storage drawer – reversible	✓	✓	✓	✓	-	-

Patient Experience	proglide300+	proglide300				
	4040-650-300	4070-650-300	4040-650-200	4070-650-200	4040-650-100	4070-650-100
Removable, Seamless Upholstered Cushion with Antimicrobial Guard	✓	✓	✓	✓	✓	✓
One touch, easy-glide stirrups with 4 lateral positions and friction lock length adjustment	✓	✓	✓	✓	✓	✓
Comfort front leg cushion	✓	✓	✓	✓	✓	✓
Heavy gauge steel extension leg pull-out to accommodate all patient heights	✓	✓	✓	✓	✓	✓

Electrical	proglide300+	proglide300				
	4040-650-300	4070-650-300	4040-650-200	4070-650-200	4040-650-100	4070-650-100
Hospital grade receptacle 110V with child safety cover (patient left)	✓	✓	✓	✓	-	-
Drawer warmer 110V	✓	✓	-	-	-	-

Quality	proglide300+	proglide300				
	4040-650-300	4070-650-300	4040-650-200	4070-650-200	4040-650-100	4070-650-100
3 year limited warranty	✓	✓	✓	✓	✓	✓
Heavy steel construction	✓	✓	✓	✓	✓	✓
Easy-clean powder coat enamel surfaces	✓	✓	✓	✓	✓	✓