



ForceTriad™ Energy Platform

KEY FEATURES

- TissueFect™ sensing technology is a Valleylab™ control system designed to precisely manage energy delivery, creating a range of options for desired tissue effect
- Enhanced LigaSure[™] tissue fusion technology permanently fuses vessels up to and including 7 mm, lymphatics, pulmonary vasculature and tissue bundles
- Automatic instrument recognition
- All-in-one unit is compatible with conventional electrosurgical instruments and all current and new LigaSure[™] instruments
- Fusion cycles are faster than the original LigaSure[™] vessel sealing generator
- Valleylab[™] mode provides a unique combination of monopolar hemostasis and dissection

BENEFITS

- Easy to use touch screen setup
- Reduces demands on staff during procedures
- ForceTriad[™] energy platform provides the surgeon improved tissue effect with consistent power delivery
- Valleylab[™] mode provides equal or superior hemostasis versus standard Cut or Blend
- Surgeons can potentially shorten their O.R. procedure time by having control of power and mode adjustments from the sterile field

OVERVIEW

Covidien, the world leader in electrosurgery for more than forty years, has integrated its latest advancements in the next generation electrosurgical device.

The Force TriVerse™ electrosurgical device allows you to access the Valleylab™ mode for a unique combination of monopolar hemostasis and dissection while using a lower power setting resulting in less char, less thermal spread and less arcing than a traditional Coagulation mode.

It is designed with you, the surgeon, in mind and offers a range of options at your fingertips for precise control and tissue effect.



ORDER INFORMATION

ForceTriad	ForceTriad™ energy platform	1 each
FT900	ForceTriad™ cart	1 each
FT950	ForceTriad [™] cart front panel	1 each
FT990	ForceTriad™ cart drawer	1 each

For more information about Valleylab™ and LigaSure™ products, visit: www.covidien.com and search for these products under "Our Brands"

To schedule an evaluation of these products, call your Covidien Sales Professional at: 800.255.8522

Covidien is a leading global healthcare products company that creates innovative medical solutions for better patient outcomes and delivers value through clinical leadership and excellence. Covidien manufactures, distributes and services a diverse range of industry-leading product lines, including Surgical Devices, Energy-based Devices, Respiratory and Monitoring Solutions, Patient Care and Safety Products, Imaging Solutions, Pharmaceutical Products and Medical Supplies. With 2006 revenue of nearly \$10 billion, Covidien has more than 43,000 employees worldwide in 57 countries, and its products are sold in over 130 countries. Please visit www.covidien.com to learn more about our business.



COVIDIEN, COVIDIEN with Logo and ™ marked brands are trademarks of Covidien AG or its affiliate. © 2007 Covidien AG or its affiliate. All rights reserved.

reorder no: LFTSYSTEMPSS 1008840 Rev. 2007/12